



AN ACCURATE AND ROBUST HANDHELD NIR



- ✓ PROTEIN,
- ✓ ASH,
- ✓ MOISTURE
- ✓ STARCH DAMAGE,
- ✓ SEDIMENTATION VALUE,
- ✓ WATER ABSORPTION,
- ✓ GLUTEN,

NIR2GO is a Diode Array (MEMS) NIR spectrometer which can be carried with you wherever you go, why waste time taking samples to a laboratory when you can have the laboratory with you.

Characteristics of the product:

Weight 190g

Battery Life Powered by tablet/laptop

Bulb Life 50,000 hrs (5+years use)

DESIGNED TO GO WHERE NO OTHER NIR INSTRUMENT HAS YET DARED TO GO



Upload your results via the internet to a server to allow instant access to the values by key staff involved the manufacturing. This means that changes can be made quickly if quality issues are found and therefore reduce the possibility of waste or even product rejection.



With support from Branscan more calibrations can be added such as Fat, Oil, Sugar and Cellulose, or if you prefer you can generate your own, this device really does remove the restrictions set in the laboratory.

A robust design with a simple interface that runs on Computer / Laptop / Tablet that is running Microsoft Windows 7/10. Really does mean that this device can be used by anybody, anywhere.



Unit 15, Lowesmoor Wharf, Worcester,
Worcestershire, WR1 2RS - England
TEL +44 (0)190529367
www.branscan.com
sales@branscan.com